

Horizontal electrophoresis apparatus

ECOSUB

Regular model

(Cat. No. E2200071)

with UV-opaque tray

(Cat. No. E2200072)

with UV-transparent tray

Regular model has standard electrodes and slim connectors

ECOSUB-D

Deluxe model

(Cat. No. E2200073)

with UV-opaque tray

(Cat. No. E2200074)

with UV-transparent tray

Deluxe model has thicker electrodes and stout connectors

- *Uni-directional migration*
- *For DNA and protein separations*
- *Gel size: 7 (w) x 10 (l) cm*
- *No. of wells: 8 x (2 combs)*
- *Buffer requirement: 200 mL*

Ordering description

ECOSUB: Horizontal apparatus with regular thickness electrodes, 2 mm male connectors, uni-directional migration & **UV-opaque tray**, gel size 7 (w) x 10 (l) cm.

ECOSUB: Horizontal apparatus with regular thickness electrodes, 2 mm male connectors, uni-directional migration & **UV-transparent tray**, gel size 7 (w) x 10 (l) cm.

ECOSUB-D: Horizontal apparatus with 1.7 times thicker electrodes, 4 mm male connectors, uni-directional migration. & **UV-opaque tray**, gel size 7 (w) x 10 (l) cm.

ECOSUB-D: Horizontal apparatus with 1.7 times thicker electrodes, 4 mm male connectors, uni-directional migration. & **UV-transparent tray**, gel size 7 (w) x 10 (l) cm.

Ready-to-use system includes

1. Moulded tank fitted with removable electrodes and male connectors
2. UV-opaque/UV-transparent tray for gel size 7 (w) x 10 (l) cm with slots for positioning two combs
3. Two PMMA combs (8 wells, 1.5 mm thick)
4. Lid
5. Electrical cords with male & female connectors designed for user safety (1 pair)
6. Instruction manual

Suggestion

Order along with Ready-to-use system an extra trays (E2203101, E2203102, E2103103, E2103104), extra combs (E2210082, E2110042), Gel slicer (E6305150)



Refer to page 12 for full list of spares and accessories