
Horizontal electrophoresis apparatus

TUBBYSUB

Regular model

(Cat. No. E2300151)

with UV-opaque tray

(Cat. No. E2300152)

with UV-transparent tray

Regular model has standard electrodes and slim connectors

TUBBYSUB-D

Deluxe model

(Cat. No. E2300153)

with UV-opaque tray

(Cat. No. E2300154)

with UV-transparent tray

Deluxe model has thicker electrodes and stout connectors

- **Uni-directional migration**
- **For DNA and protein separations**
- **Gel size: 15 (w) x 10 (l) cm**
- **No. of wells: 15/20 x (2 combs)**
- **Buffer requirement: 500 mL**

Ordering description

Tubbsub: Horizontal apparatus with regular thickness electrodes, 2 mm male connectors, uni-directional migration & **UV-opaque** tray, gel size 15 (w) x 10 (l) cm.

Tubbsub: Horizontal apparatus with regular thickness electrodes, 2 mm male connectors, uni-directional migration & **UV-transparent** tray, gel size 15 (w) x 10 (l) cm.

Tubbsub-D: Horizontal apparatus with 1.7 times thicker electrodes, 4 mm male connectors, uni-directional migration. & **UV-opaque** tray, gel size 15 (w) x 10 (l) cm.

Tubbsub-D: Horizontal apparatus with 1.7 times thicker electrodes, 4 mm male connectors, uni-directional migration. & **UV-transparent** tray, gel size 15 (w) x 10 (l) cm.

Ready-to-use system includes

1. Tank fitted with removable electrodes and male connectors
2. UV-opaque tray/UV-transparent tray for gel size 15 (w) x 10 (l) cm with slots for positioning two combs
3. Two PMMA comb (15 wells, 1.5 mm thick)
4. Lid
5. Electrical cords with male & female connectors designed for user safety (1 pair)
6. Instruction manual



Suggestion

Order along with Ready-to-use system: extra trays (E2103104, E2203102, E2303106), extra combs (E2110042, E2210082, E2310152, E2310202)

Refer to page 12 for full list of spares and accessories
